

Tween[®] 20 (Polysorbate 20)

Tween[®] 20 is a polyoxyethylene sorbitol ester that belongs to the polysorbate family. It is a nonionic detergent widely used in biochemical applications.

Request quote



Product Code	Product Description	Pack Size	Cas No.	HSN Code	List Price (INR)
Q28625	Tween 20 (Polysorbate 20)	500 ml	9005-64-5	34021300	1666

TWEEN 20 is Used for solubilizing membrane proteins during isolation and purification. Polysorbate 20 is used as a wetting agent in flavored mouth drops. It acts as a washing agent in Western blots and ELISAs and prevents antibody binding. It is a solubilizing agent of membrane proteins. In pharmaceutical chemistry, it is used as an excipient to stabilize emulsions and suspensions. Further, it is used as a dispersing agent, stabilizer, lubricator and also used in creams, salves, ointments, protective creams, balms, pomades, lipsticks, mascaras and glosses.

